

PAJUNK®

MonoBalloon Systems
Fully Disposable Balloons for
Minimally Invasive Procedures




MADE IN GERMANY

Dissecting balloons


The Pajunk dissecting balloon systems allow for optimum, fast and effective dissection of the preperitoneal space with direct visualisation.

The offering includes an oval balloon for bilateral hernias and round balloon for unilateral hernia repairs.

Unilateral hernia repairs



**Bilateral
hernia repairs**



Green marking
→ Safe application
proceeding

Obturator

- Low profile
- Blunt tip



Access balloons


Pajunk offers two different balloon systems to secure the access:
ring-anchor and structural balloon systems.

MonoBalloon have a tried and tested standard seal system with a silicone valve,
and VarioBalloon systems have a universal double seal system.

MonoBalloon system

with standard seal

Systems are available in short and long lengths (e.g. for LERP).



MonoBalloon SL

- ➔ Easier insertion
- ➔ Higher stability

➔ All balloons are latex-free

VarioBalloon system

with universal double seal

- Completely gas-tight
- Secure bayonet lock
- Easy twist lock


Obturator Mono

Mono PROTECT

Obturator Vario

Vario PROTECT

Ø3 – Ø5.8mm
Ø5 – Ø11mm



NEW

- Check valve**
- Automatic system closure
 - One-hand operation
 - Quick & easy balloon filling

VarioBalloon SL

VarioBalloon URO

VarioBalloon RA

→ Trauma reduction with blunt tip

→ All balloons are latex-free

→ Trauma reduction with blunt tip

MonoBalloon systems



Dissecting balloons

Product	Size	Item no.	PU	Ball pump
MonoBalloon DIL-XL Unilateral dilation balloon (round, extra large), sterile	Ø 11 mm	1284-60-11	5	•
MonoBalloon DIL-BI Bilateral dilation balloon (oval), sterile	Ø 11 mm	1284-50-11	5	•



Access balloons

Product	Size	Item no.	PU	Ball pump	30 ml syringe
MonoBalloon SB Structural balloon and instrument channel, sterile	Ø 11 mm	1284-20-11S	5	•	
MonoBalloon RA Ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1284-30-11S	5	•	
	Ø 5.8 mm	1284-30-58S	5	•	
MonoBalloon RA PROTECT Ring-anchor balloon and instrument channel, sterile	Ø 5.8 mm	1284-31-58S	5	•	
MonoBalloon RA-S Short, ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1284-40-11S	5	•	
MonoBalloon URO Flat ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1284-70-11S	5	•	
	Ø 5.8 mm	1284-70-58S	5	•	
MonoBalloon URO PROTECT Flat ring-anchor balloon and instrument channel, sterile	Ø 5.8 mm	1284-71-58S	5	•	
MonoBalloon URO-S Short, flat ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1284-75-11S	5	•	
MonoBalloon SL Flat ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1284-80-11S	5	•	
MonoBalloon SL-S Short, flat ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1284-85-11S	5	•	

MonoBalloon systems – sets

Double balloon sets



MonoBalloonDuo DIL-XL + RA

Unilateral dilation balloon (round, extra large), ring-anchor balloon and instrument channel, sterile

Size

Item no.

PU

Ball pump
30 ml syringe


Ø 11 mm

1284-60-20S

5

•

•



MonoBalloonDuo DIL-XL + SB

Unilateral dilation balloon (round, extra large), structural balloon and instrument channel, sterile

Ø 11 mm

1284-60-24S

5

•



MonoBalloonDuo DIL-XL + URO

Unilateral dilation balloon (round, extra large), URO ring-anchor balloon and instrument channel, sterile

Ø 11 mm

1284-60-28S

5

•

•



MonoBalloonDuo DIL-BI + RA

Bilateral dilation balloon (oval), ring-anchor balloon and instrument channel, sterile

Ø 11 mm

1284-50-20S

5

•

•



MonoBalloonDuo DIL-BI + SB

Bilateral dilation balloon (oval), structural balloon and instrument channel, sterile

Ø 11 mm

1284-50-24S

5

•



MonoBalloonDuo DIL-BI + URO

Bilateral dilation balloon (oval), URO ring-anchor balloon and instrument channel, sterile

Ø 11 mm

1284-50-28S

5

•


•

VarioBalloon systems with universal double seal



Access balloons

Product	Size	Item no.	PU	30 ml syringe	20 ml syringe	10 ml syringe
VarioBalloon RA						
Ring-anchor balloon, instrument channel and obturator, sterile	Ø 11 mm	1283-30-11S	5	•		
	Ø 5.8 mm	1283-30-58S	5			•
VarioBalloon RA PROTECT						
Ring-anchor balloon, instrument channel and PROTECT-obturator, sterile	Ø 5.8 mm	1283-31-58S	5			•
VarioBalloon URO						
Flat ring-anchor balloon, instrument channel and obturator, sterile	Ø 11 mm	1283-70-11S	5	•		
	Ø 5.8 mm	1283-70-58S	5		•	
VarioBalloon URO PROTECT						
Flat ring-anchor balloon, instrument channel and PROTECT-obturator, sterile	Ø 5.8 mm	1283-71-58S	5		•	
VarioBalloon SL						
Flat ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1283-80-11S	5		•	
VarioBalloon SL-S						
Short, flat ring-anchor balloon and instrument channel, sterile	Ø 11 mm	1283-85-11S	5		•	



PAJUNK® GmbH
 Medizintechnologie
 Karl-Hall-Strasse 1
 D-78187 Geisingen/Germany
 Phone +49 (0) 77 04/92 91-0
 Telefax +49 (0) 77 04/92 91-6 00
www.pajunk.com

PAJUNK® Medical Produkte GmbH
 D.A.CH • BeNeLux
 Karl-Hall-Strasse 1
 D-78187 Geisingen/Germany
 Phone +49 (0) 77 04/80 08-0
 Telefax +49 (0) 77 04/80 08-150
www.pajunk.com