

MultiStim ECO

Nerve Stimulator
for Nerve Blocks



MultiStim ECO

Dual Guidance Nerve Stimulator for Nerve Blocks

The MultiStim ECO nerve stimulator was developed primarily as a Dual Guidance stimulation device, however, it can also be used for stimulation-supported puncture techniques. The easy-to-use, compact device has been optimised by integrating the most frequently used functions and parameters. This makes handling fast, simple, and more cost-effective. The needle can be connected to the device with no additional patient cables required.

- Limited, easy settings
- Snap-on button
- Direct connection to the needle



Fixation Directly on the Patient

- Adhesive electrode underneath the device with a snap-on button connection
- Increased convenience



Direct Needle Connection

- No patient lead necessary
- Increased convenience



Variation of the Current Strength

- Six-level amplitude scale from 0.2 to 2.0 mA
- Arrow buttons are used for the plus/minus key



Fixed Defined Settings

- Frequency: 1 Hz/
stimulation pulsewidth: 0.1 ms



Adaptable Patient Cable (PlugX)

- Allows for fixation of the ECO away from the patient
- Extended cable length

MultiStim ECO	Description	Item no.	PU
Stimulator	with device case	1151-94-50	1
	PlugX – patient cable with plug 2.0mm for needle and clamp adapter for electrode	1151-94-21	1
Accessories	Extension cable for connecting stimulation needles, autoclavable	01151-861F	10
	Extension cable for connecting stimulation needles and stimulators, disposable	01151-861Q	10

PAJUNK® GmbH Medizintechnologie
Global Headquarters
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

PAJUNK® Medical Produkte GmbH
Central Europe
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de & pajunk.eu

PAJUNK® Swiss AG
Switzerland
Neuhofstrasse 5A
6340 Baar, Switzerland
Tel.: +41 417 609 627
info@pajunk.ch
pajunk.ch

PAJUNK® Medical System L.P.
United States of America
4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 888 972 5865
Fax: +1 678 514 3388
info@pajunk-usa.com
pajunkusa.com

PAJUNK® UK Medical Products Ltd
United Kingdom
Unit C1 The Waterfront
Goldcrest Way, Newburn Riverside
Newcastle upon Tyne, NE15 8 NY, UK
Tel.: +44 191 264 7333
info@pajunk.co.uk
pajunk.co.uk