


SonoBlock II

*Echogenic, Non-Stimulating
Single Shot Plane Block Needle*

SonoBlock needles were developed by PAJUNK® for single shot nerve block applications without stimulating, allowing for placement under ultrasound guidance only. The needles also includes the patented Cornerstone Reflectors for optimized needle echogenicity.



¹ Edgcombe H, Hocking G. (2010 March–April) Sonographic identification of needle tip by specialists and novices: a blinded comparison of 5 regional block needles in fresh human cadavers. Reg. Anesth. Pain Med., 35(2): 207–211.

SonoBlock II

Echogenic, Non-Stimulating
Single Shot Plane Block Needle

SonoBlock II	Size	Item no.	NRFit® Item no.	PU
Facet tip	22G x 50mm (2")	001280-74	001260-74	10
	22G x 80mm (3 1/5")	001280-71	001260-71	10
	21G x 100mm (4")	001280-77	001260-77	10

SonoBlock II	Size	Item no.	NRFit Item no.	PU
Facet S tip	22G x 50mm (2")	001281-74		10
	22G x 80mm (3 1/5")	001281-71		10
	21G x 100mm (4")	001281-77		10
	20G x 120mm (4 3/4")	001281-72		10

NRFit® is a registered trademark of GEDSA used with their permission.

PAJUNK® GmbH Medizintechnologie
International Sales
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

PAJUNK® Medical Produkte GmbH
Central Europe
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de & pajunk.eu

PAJUNK® Swiss AG
Switzerland
Neuhofstrasse 5A
6340 Baar, Switzerland
Tel.: +41 417 609 627
info@pajunk.ch
pajunk.ch

PAJUNK® Medical System L.P.
United States of America
4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 888 972 5865
Fax: +1 678 514 3388
info@pajunk-usa.com
pajunkusa.com

PAJUNK® UK Medical Products Ltd
United Kingdom
Unit D1 The Waterfront
Goldcrest Way, Newburn Riverside
Newcastle upon Tyne, NE15 8NY, UK
Tel.: +44 191 264 7333
info@pajunk.co.uk
pajunk.co.uk